



Thermo Scientific iCAP RQ ICP-MS



Simplicity, productivity and robustness
for routine labs

Thermo
SCIENTIFIC



Thermo Scientific iCAP RQ ICP-MS

The reliability, analytical performance and ease of use needed to meet the demands of the highest throughput labs

More productivity

Choose the Thermo Scientific™ iCAP™ RQ ICP-MS for a complete multi-element analysis solution for your high-throughput routine lab. With comprehensive interference removal for assured data accuracy and intuitive workflows to boost productivity, this single quadrupole (SQ) ICP-MS will expand your analytical capabilities.

Our user-inspired hardware and software maximizes ease-of-use and streamlines workflows to achieve 'right-first-time' analysis, essential to all busy laboratories.

Increase throughput, minimize sample re-runs and save analysts' valuable time and effort with intuitive features and broad compatibility with automated accessories.

- 'Get Ready' one-click set-up for easy, reproducible analysis
- Cone change in as little as 60 seconds
- Fast, easy operation with quick-connect sample introduction components

More robustness

With its logical, problem-solving hardware design that meets real laboratory needs through minimized maintenance and maximized up-time, the iCAP RQ ICP-MS is a powerful, reliable work-horse analyzer.

Run samples 24/7 with complete confidence.

- Matrix tolerance and minimized drift in challenging samples, such as seawater
- Reliable, robust plasma operation, even for volatile organic solvents such as 100% acetonitrile
- Advanced heat management for long-term mass stability
- All new electronics and RF generator design



Bench height, open-access, sample introduction system.



Integrated HD camera captures plasma images during all stages of operation.






Why Choose SQ-ICP-MS for Your Routine Laboratory?


For the simplest possible method development and smart design features that drive productivity

- Measure all analytes in a single helium Kinetic Energy Discrimination (He KED) mode
- In He KED mode, even lithium and beryllium can be measured at sub ppt levels
- Excellent interference removal with high ion transmission
- Reactive gas capability targets specific interference ions for challenging applications

Interface	When	Why	Who
Robust	Long-term analysis of dirty, high matrix samples	Minimal drift and reduced user maintenance	Environmental and routine industrial applications
High Matrix	Typical, everyday matrix types	Best balance for sensitivity and matrix tolerance	Food safety, pharmaceutical and drinking water analysis
High Sensitivity	Advanced applications; e.g. laser, nanoparticles	Best signal to noise ratio and lowest detection	Geoscience, advanced materials, research




Choose from three dedicated interface types designed for maximum coverage of sensitivity and dynamic range.




Quadrupole set to filter out exact mass of target analyte.

QCell flatapole dynamically applies LMCO relative to target analyte.



QCell low mass cut off (LMCO) removes precursors to interferences and reduces BECs even further than He KED alone.



Compact, vertical design to save valuable lab space.

Powerful detector for high sensitivity.

QCell combined with unique flatapole design enables high data accuracy.

Simplified methods with He KED comprehensive interference removal.

Zero user maintenance beyond the slide valve.

New solid state RF generator and all new electronics ensures stability and maximizes uptime.

Plasma robustness to cope with the most challenging matrices.

Quick connect, self-aligning sample introduction system.

Drop-down door for simplified cone access.



Qtegra Intelligent Scientific Data Solution (ISDS) Software

Universal software for ultimate ease-of-use and productivity in high -throughput labs and demanding applications

Maximize productivity and simplify workflows utilizing intelligent features and automation.

Method development is streamlined with Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution™ (ISDS) software's intuitive user features; productivity is enhanced through a range of automated sample handling options, all managed within a logical dashboard interface.

Sample Handling

Our analysis-specific sample handling kits and additional mass flow controllers simplify instrument configuration and enable reproducibility.

Plug-in compatibility with a wide range of industry standard autosamplers and sampling valves optimizes your workflows.


Auto-dilution

Increase throughput and automate the entire workflow from sample preparation to data review with an intelligent sample management accessory. Ensure error free operation of the full sample batch, including:

- Complete calibration from a single standard
- Intelligent dilution of high matrix or over range samples

Argon Gas Dilution

Eliminate the need for manual sample dilution and reduce matrix build-up on the interface by using automated Argon Gas Dilution (AGD). This feature provides a simple solution for direct analysis of high matrix samples at % level concentrations, enabling extended periods of operation without cleaning.



Accelerate your innovation with the advanced applications enabled by **hyphenated techniques**.

Our simple coupling kits and a range of Qtegra ISDS software plug-ins allow easy pairing with IC/LC, GC and laser ablation systems.

Explore the expanding field of **nanoparticle analysis** with our npQuant plug-in for Qtegra ISDS software.

Reliable, unattended operation helps deal with sample backlogs.



Our easy-to-use Qtegra ISDS software delivers all the support features essential to any busy laboratory, while containing all the flexibility needed to achieve the most challenging applications.

This software approach, common to other Thermo Scientific ICP techniques, reduces training requirements and makes the addition of a SQ-ICP-MS into your lab, a simple process.

Low Cost per Sample

Operator time is freed up with simple workflows and a minimal number of steps per task.

The 'Get Ready' feature takes your iCAP RQ ICP-MS from standby to ready for analysis, through a fully automated process, saving time and ensuring consistent analysis.

A fully integrated QA/QC protocol reduces failures, driving down re-runs and cost per sample.

Data and Report Management

Automated report generation and export features, including full LIMS compatibility, makes data management effortless.


Customize report formats to fit with the needs of the laboratory or your customers.

Compliance Support

All the features needed to enable compliance with regulatory requirements are available within Qtegra ISDS software:

- EPA and FDA regulatory compliance and audit support
- Data security and access control
- Compliance management tools

Recovery results:				
User name	Index	Label	51V (KED)	60Ni (KED)
MP14-PC\MP14 - second user	6	QL-1	101.811 %	103.496 %
MP14-PC\MP14 - second user	7	QL-2	99.219 %	101.195 %
MP14-PC\MP14 - second user	8	QL-3	100.279 %	102.939 %
MP14-PC\MP14 - second user	9	QL-4	100.180 %	103.087 %
MP14-PC\MP14 - second user	10	QL-5	100.677 %	103.699 %
MP14-PC\MP14 - second user	11	QL-6	98.594 %	102.038 %
Average		RSD	100.4712 %	102.7407 %
			1.0544 %	0.7639 %
			99.694 %	100.407 %
			100.427 %	102.909 %
			101.375 %	103.373 %
			99.907 %	101.901 %
			97.024 %	99.432 %
			102.300 %	104.294 %
			100.1211 %	102.0534 %
			1.7997 %	1.8131 %
			Pooled Recovery	100.1461
			Pooled RSD	1.4064 %
				101.8619
				101.3727
				1.5266 %
				1.3309 %



Advanced reporting capability.

Growing need for nanoparticle characterization

Enhance food and environmental nanoparticle characterization
with the npQuant plug-in for accurate single particle
ICP-MS (spICP-MS) analysis.

In Action: Environmental, Food Safety, Pharmaceuticals, Clinical Research

Handle a broad variety of sample types with ease and accuracy

Environmental

Confidently measure trace and major analytes in some of the most challenging matrices. Optimize productivity and accuracy with comprehensive interference removal for reliable multi-elemental analysis.

Whether analyzing drinking water, waste water or soil digests, the high-throughput, iCAP RQ ICP-MS meets demanding environmental requirements, such as US EPA 200.8 with ease.



More certainty with speciation

Species information is critical in environmental, drinking water, food and pharmaceutical analysis. Easily couple your iCAP RQ ICP-MS with an IC system for assured species data using IC-ICP-MS.



Pharmaceutical Compliance

Comfortably meet the most demanding global pharmaceutical regulations and legislation, including;

- International Council on Harmonization (ICH) guideline Q3D
- U.S. Pharmacopoeia (USP) chapters, 232, 233 and 2232

Comply with the most rigorous data audit and security measures. Our Qtegra ISDS software is fully compliant with Food and Drug Administration (FDA) 21 CFR Part 11 and is equipped with full IQ/OQ procedures for simple implementation in GMP/GLP regulated environments.



Clinical Research

Detect the smallest traces of toxic and essential elements, with the robustness, sensitivity and LoDs needed for complete analyte quantitation, even in tough matrices.

Accelerate your research with assured accuracy, utilizing analysis-boosting features such as auto-optimization of QCell reaction mode.

Food Safety

Simultaneously measure toxic and essential elements for food quality and safety assessment. Be ready to meet food safety legislation and regulatory requirements with a suite of integrated QC features.

Couple metal-free Ion Chromatography systems to ICP-MS to easily and accurately speciate critical elements such as chromium, arsenic and mercury.



In Action: Materials, Nuclear, Geoscience Exploration

Robust, reliable performance, even in the toughest environments

Solid Sampling with Laser Ablation

Expand capabilities, minimize sample preparation and improve analysis times for difficult samples. Easily integrate laser ablation systems to obtain the best spatial and depth characterization.

Advantages of the iCAP RQ ICP-MS

- Analyze challenging matrices, across a broad range of sample types
- Reliable and robust plasma, even in the presence of 100% organic solvents
- Effective comprehensive interference removal, for assured accuracy
- Simple maintenance, minimal down-time
- Straightforward to set up, ergonomic and intuitive to use
- Common software platform across our ICP techniques reduces training and enables operator flexibility
- Easily configured for laser ablation solid sampling
- Ability to analyze natural and engineered nanoparticles
- Hardware and software ready to easily integrate advanced hyphenated techniques



Material Analysis

Whether running low level alloy QA/QC analysis or qualifying advanced materials in the aerospace industry, the iCAP RQ ICP-MS can easily be configured to your needs.

Simplify material analysis with fully integrated, fast sampling accessories and intelligent dilution for automated sample preparation.

Flexibility of collision/reaction modes for optimal interference removal, enables accurate analysis in a broad range of material types.



Nuclear

Achieve advanced isotopic detection utilizing a high sensitivity interface and superior mass stability. Unique ion optics offer low background and high signal to noise ratios for optimal detection limits.

Collisional focusing with our unique QCell technology ensures the abundance sensitivity required for the challenge of measuring isotope ratios.



Geoscience Exploration

Robust long-term analysis of the toughest rock digests and lithium metaborate fusions. Customized interface ensures minimal drift over time and integrated QC checks enable automatic self assessments of performance.



We manage your instruments so you can focus on the science

When you buy a Thermo Scientific product, you gain the peace of mind that comes from the backing of a global team of Unity Lab Services experts committed to your long-term success.



Thermo Scientific™ iCAP™ RQ ICP-MS, part of the iCAP™ Qnova™ Series ICP-MS

A total, multi-element solution for high-throughput routine laboratories.

Expand your laboratory portfolio with this powerful, compact analyzer. With uniquely reliable cell mode performance and user-friendly platform software, SQ-ICP-MS analysis is easier than ever. Ideal for labs with a focus on environmental, food, water and geoscience analyses.



Compact design



Easy installation



User-friendly operation

Simplicity, productivity and robustness
for routine labs

www.thermofisher.com/SQ-ICP-MS

©2016 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. ESI prepFAST 2 is used in trade by Elemental Scientific Inc. LSX 213 G2+ and ASX 560 are used in trade by Teledyne CETAC Technologies. USP is a trademark of United States Pharmacopeial Convention. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. For research use only. Not for use in diagnostic procedures.

BR43317-EN 0216C

Thermo
S C I E N T I F I C

A Thermo Fisher Scientific Brand